

# V200

lightweight

Removable wheels



In 3 heights and depths adjustable  
armpad

Foldable and  
detachable arm rests











Foldable and  
detachable  
foot supports

Adjustable seatdepth

Adjustable footplates



## Technical information V200

									
total width	seat width	height armrests	seating height	total height	total length	seating depth	height back	weight	users weight
570	390	220-230-240	470-500-530	850-880-910-940	1010	440-460	400	15,2	130
600	420	220-230-240	470-500-530	850-880-910-940	1010	440-460	400	heaviest part 9	
620	440	220-230-240	470-500-530	850-880-910-940	1010	440-460	400		
640	460	220-230-240	470-500-530	850-880-910-940	1010	440-460	400		
660	480	220-230-240	470-500-530	850-880-910-940	1010	440-460	400		
680	500	220-230-240	470-500-530	850-880-910-940	1010	440-460	400		
710	530	220-230-240	470-500-530	850-880-910-940	1010	490	400		
740	560	220-230-240	470-500-530	850-880-910-940	1010	490	400		
780	600	220-230-240	470-500-530	850-880-910-940	1010	490	400		





Complies to crash test ISO 7176-19

Sizes in mm, weight in kg. Deviations in sizes and weight (1,5 mm, 1,5 kg, 1,5°) are possible.

### Model V200

The V200 is a light weight manually propelled wheelchair for temporary or permanent use with optimised settings.

- 4 seat heights
- 2 seat depths
- 3 armrests heights
- 3 armpad depths
- 5 seat inclinations 0° - 10°
- a range of available accessories

Frame	Upholstery
	
Champagne C59	Black Nylon/Skai
	
Anodic Brown C60	Black Nylon/Skai

### The wheelchair must be adapted to you, not vice versa.

Our basic principle is: the wheelchair must be adapted to you, not vice versa. This is why our wheelchairs have a modular design. All types of wheelchairs can be adapted, have numerous accessories and all options are mutually interchangeable. This is the result of our continuous research and more than 53 years of experience.

### Your guarantee of quality.

Regarding quality, the **V200** complies with the highest norms, laid down by certified test institutes in most European countries.

## N.V. Vermeiren N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

Tel.: +32 (0) 3 620 20 20  
Fax: +32 (0) 3 666 48 94

[www.vermeiren.com](http://www.vermeiren.com) - [info@vermeiren.be](mailto:info@vermeiren.be)



**VERMEIREN**  
Opening up new horizons in mobility

For more information, contact your local dealer

