

201 Delta

Rollator



201 Delta is a foldable 3-wheel rollator suitable for indoor as well as outdoor use. With the included basket and storage bag, it can easily be used while shopping. Furthermore, this rollator is equipped with height adjustable handles with brakes. Because of its foldable lightweight frame, the **201 Delta** rollator is easy to transport.

total width	680 mm
total height	820 - 920 mm
total length	576 mm
total weight	5,8 kg
max. users weight	125 kg
wheels (front/rear)	8 x 1 ^{1/4} ø 190 mm
colour(s)	champagne

Measurement tolerances: ±15mm / 1,5kg.
We reserve the right to introduce technical changes.



N.V. VERMEIREN N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium
Tel.: +32 (0)3 620 20 20 - Fax: +32 (0)3 666 48 94 - info@vermeiren.be - website: www.vermeiren.com

Model shown only for illustration purposes. Certain features may be optional.

Opening up new horizons in mobility