



# V300

## Lightweight Active Wheelchair

The V300 is a lightweight manual wheelchair for temporary or permanent use with a fixed backrest. What makes the V300 unique is its seat depth of up to 52 cm. The largest seat depth in our range, providing maximum comfort to its user. The adaptable hook and loop backrest and seat cushion with comfortable padding are always included. V300 also has other options available such as a headrest, table and side supports.



**VERMEIREN**  
we care for you

## Standard

- Available in 9 different sizes (from 39 cm to 60 cm\*)
  - \* Sizes 53 - 56 - 60 cm: XL version
- 2 colour options (champagne and blue)
- Tension adjustable backrest and seat
- Armrests adjustable in height and depth (3 x 2 adjustments)
- Removable swing-away leg rests
- Removable armrests for easier transfers

## Colours



## Features



*Armrests adjustable in 3 different heights and 3 different depth settings*



*Backrest and seat are adjustable in tension*



*Adjustable seat depth*



*Numerous possibilities of adjustments and configurations*

## Technical data

■ Total width	600/630/650/670/690/710 mm
■ Total height	850 - 940 mm
■ Total length	1040 mm
■ Folded width	310 mm
■ Seat width	390/420/440/460/480/500 mm
■ Height armrests	220 - 240 mm

■ Seat height	440 - 530 mm
■ Height backrest	420 mm
■ Seat depth	440 - 520 mm
■ Weight	17,7 kg
■ Maximum user weight	130 kg
■ Additional info	

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **V300** complies with the highest norms, laid down by certified test institutes in most European countries.

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



**For more information, contact your local dealer**