

# Vera

## Transfer board



**Vera is a plastic transfer board with a 360° rotating disc.** Equipped with multiple holes to allow the water to drain and prevent slippage. No tools are required and easy to install on most bathtubs available on the market.

■ Total width	350 mm
■ Total length	730 mm
■ Total height	175 mm
■ Diameter swivel disk	325 mm
■ Weight	3,75 kg
■ Max. user weight	130 kg
■ Color(s)	white
■ Additional information	CE

Measurement tolerances:  $\pm 15$ mm / 1,5kg.  
We reserve the right to introduce technical changes. Models shown only for illustration purposes. Certain features may be optional.

