

V100

lightweight

Removable wheels

In 3 heights and depths adjustable
armpad











Foldable and
detachable arm rests

Foldable and
detachable
foot supports

Adjustable seatdepth



Technical information V100

									
570	390	220-230-240	470-500-530	880-910-940	1010	430-450	420	16,7	130
600	420	220-230-240	470-500-530	880-910-940	1010	430-450	420	heaviest part 10,5	
620	440	220-230-240	470-500-530	880-910-940	1010	430-450	420		
640	460	220-230-240	470-500-530	880-910-940	1010	430-450	420		
660	480	220-230-240	470-500-530	880-910-940	1010	430-450	420		
680	500	220-230-240	470-500-530	880-910-940	1010	430-450	420		
710	530	220-230-240	470-500-530	880-910-940	1010	490	420		
740	560	220-230-240	470-500-530	880-910-940	1010	490	420		
780	600	220-230-240	470-500-530	880-910-940	1010	490	420		





Complies to crash test ISO 7176-19

Sizes in mm, weight in kg. Deviations in sizes and weight (1,5 mm, 1,5 kg, 1,5°) are possible.

Model V100

The V100 is a light weight manually propelled wheelchair for temporary or permanent use with a range of settings.

- 3 seat heights
- 2 seat depths
- 3 armrests heights
- 3 armpad depths
- a range of available accessories

Frame	Upholstery
	
Metallic Blue C58	Black Nylon/Skai
	
Metallic Inox C55	Black Nylon/Skai

The wheelchair must be adapted to you, not vice versa.

Our basic principle is: the wheelchair must be adapted to you, not vice versa. This is why our wheelchairs have a modular design. All types of wheelchairs can be adapted, have numerous accessories and all options are mutually interchangeable. This is the result of our continuous research and more than 53 years of experience.

Your guarantee of quality.

Regarding quality, the **V100** complies with the highest norms, laid down by certified test institutes in most European countries.

N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be



VERMEIREN
Opening up new horizons in mobility

For more information, contact your local dealer

