





Isometric view




Side view




Detail B



Top view



Detail C



Detail A

No.	Title	Art.nr.	Sticker	Number
1	Sticker Vermeiren 85x15	T001051		2
2	Sticker taxi hooks	T001076		4
3	Sticker serial number	T001097		1
4	Sticker name Eclips X4	T001038		2

STICKER CONFIGURATION ECLIPS X4

Plan Number

E000686/C

Scale: 3:40	Drafted	Verified	Confirmed
Date	8/02/2011	27/01/2015	27/01/2015
Initials	TT	TT	PJ

