



No.	Plan	Article	Qty
1	Y011987	1669854	1
2	Y011148	1669754	6
3	Tube cap 15x15	1005899	1
4	Lever M8x12	1900712	3
5	XP01548	1001524	3
6	Washer M8	1996214	3
7	E00066472	2998017	1
8	E00052508	2998018	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

2998010

Plan No.
E00066356/A.1

PLEGO L58			
Scale: 1:4	Drafted	Verified	Confirmed
Date	17.08.2023	17.08.2023	17.08.2023
Initials	BK	WM	TK