

Material					
No.	Article	S	M	L	XL
1	2991014	0,062	0,068	0,082	0,094
2	1009517	4,8	4,8	4,8	4,8
3	1009512	12	12	12	12
4	1009508	7	7	7	7
5	2991330	5	5	5	5
6	2991886	0,9	1,2	1,5	1,8

	Sizes [cm]		Article	
	LENGHT	FR	BK	
S	31	2991259	2991255	
M	34	2991260	2991256	
L	41	2991261	2991257	
XL	47	2991262	2991258	

Mounting			
No.	Plan	Article	Quantity
1	XP01976	2993641	4
2	XP01980	2993813	4
3	XP01979	2993812	4
4	M6X30 collar bolt	2994294	4
5	YP07178	2760099	4



**VERMEIREN
NV**

Vermeirenplein 1-15
2920 Kalmthout
www.vermeiren.com



Made in Poland

NeoFlex

**U79Neo.
size M**

contains:
90% Neopren
8% Poliester
2% Elastan

**SGS INSTITUT
FRESENIUS**



EP00346 / B

NEO U79 H-TYPE SHOULDER BELTS				Plan No.
Scale: N/A	Drafted	Verified	Confirmed	EP00346/ B
Date	31/07/2018	04/03/2019	04/03/2019	
Initials	FD	AS	TK	