



Isometric view  
Scale: 1:10

Material			
No.	Title	Article	Sticker
1	Sticker Ø35 logo Vermeiren	T001002	
2	Sticker taxi hooks	T001076	
3	Sticker serial number	T001025	
4	Sticker recyclage electro	1997820	
5	Sticker name Vermeiren new large	T001052	

SPRINGER KIDS_FRAME WITH STICKERS				Plan Number
				EP00093/A
Scale: 1:1	Drafted	Verified	Confirmed	
Date	17/09/2012	17/09/2012	17/09/2012	
Initials	AS	MK	TK	